Sanctband

Sanctband *Active*

LATEX PRODUCTS SAFETY MANUAL

Safety Instructions – Please Note

1. Natural product with reliable linear resistance behaviour. Sanctband is a functional piece of exercise device for physical therapy, muscle increase and fitness.
2. **ATTENTION:** Sanctbands are made from natural latex which can cause allergic reactions.
3. When using Sanctband, always follow the recommendations of your doctor or therapist.
4. Children may use Sanctband only under adult supervision and after a therapy-accompanying introduction. Not suitable for infants.
5. Check Sanctband for damages (holes, cracks, cuts) prior to every use. Do not use damaged bands, replace them instead.
6. Do not overstretch Sanctband ($F \geq 200\%$), do not create too much tension. Choose the band that matches your physical condition and muscle strength.
7. Be careful with the exercises at head height! Keep away from eyes, the band could slip out of hand or holder and snap in the face!
8. Ensure the band cannot release from door frame or loosen from the chair or other objects.
9. Ensure that door is shut while exercising.
10. Sanctband's Recommended Cleaning Method: Clean with wet clean cloth.

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Wujiang Ltd
Sanctband Health Care
www.sanctband.com

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