EXERCISE  
**LOOP BAND**  
INSTRUCTION MANUAL

Sanctuary Health Sdn Bhd (583967-A)  
16, Persiaran Perindustrian Kanthan 1,  
Kawasan Perindustrian Kanthan,  
31200 Chemor, Serak-Bandar Ritzuan, Malaysia.  
Tel: +605-201 2800 Fax: +605-201 7800  
Website: www.sanctuaryhealth.com.my

CE EU REP conMed GmbH  
Friedenstraße 9  
11055 Berlin  
GERMANY conMed GmbH  
Friedenstraße 9  
11055 Berlin  
GERMANY LATEX MD

Date of issue: 2025-09-08

**Safety Instructions – Please Note**

- Natural product with reliable linear resistance behaviour, Sanctband | Sanctband Active is a functional piece of exercise device for physical therapy, muscle increase and fitness.
- ATTENTION: Sanctband | Sanctband Active are made from natural latex which can cause allergic reactions.**
- When using Sanctband | Sanctband Active, always follow the recommendations of your doctor or therapist.
- Children may use Sanctband | Sanctband Active only under adult supervision and after a therapy-accompanying introduction. Not suitable for infants.
- Check Sanctband | Sanctband Active for damages (holes, cracks, cuts) prior to every use. Do not use damaged loop band, replace them instead.
- Do not overstretch Sanctband | Sanctband Active (F≥200%), do not create too much tension. Choose the loop band that matches your physical condition and muscle strength.
- Be careful with the exercises at head height! Keep away from eyes, the loop band could slip out of hand or holder and snap in the face!
- Ensure the loop band cannot release from door frame or loosen from the chair or other objects. Ensure that door is shut while exercising.
- Sanctband's | Sanctband Active's Recommended Cleaning Method: Clean with wet clean cloth.
- Shelf Life: 5 years from the date of manufacture if unopened and unused.
- Disposal of the devices:** Follow local disposal procedure as per user facility.

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